Scissors
剪刀

3.01-12

Basic Informations
基本情報

3.02 - 3.04

Scissors flat handle
剪刀(平柄)

3.05 - 3.06

Scissors round handle
剪刀(丸柄)

3.07 - 3.08

Special Scissors round handle, by Blondeel, curved
特殊剪刀(丸柄)

3.09 - 3.10

Special Scissors flat handle, angulated 45°
特殊剪刀(平柄)、角度45°

3.11

Special Scissors flat handle / round handle
特殊剪刀(平柄) / 特殊剪刀(丸柄)

3.12

S&T®

Basic informations

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Scissors

All **S&T®** microsurgical scissors are skillfully handcrafted to microsurgical scissors feature: exacting standards. You will immediately notice the advantages of these scissors:

- Flawless design due to years of experience in the field
- A uniquely **smooth cutting action** due to precise engineering
- **Minimum closing pressure** in operation due to equal pressure distribution along the blades

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S&T® Microsurgical Scissors Patterns

Dissection

This pattern features **curved** blades with tips that are rounded to avoid possible damage to structures during dissection.

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Adventitia

This pattern features **straight** blades with tips that are pointed for fine trimming and suture cutting.

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Nerve/Tendon

This pattern features straight blades with one **finely serrated** blade, ensuring a controlled even cut when preparing tissue such as nerves or tendons. Due to design characteristics, **S&T®** serrated scissors provide exceptionally smooth cutting action compared to other serrated scissors.

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Basic informations

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Dissecting or Adventitia Pattern?

Microsurgical scissors in the S&T® product offering are available in either dissecting or adventitia patterns. The distinction lies in the sharpness of the tips and the curvature of the blades. Dissection scissors are rounded at the tip with a radius of 0.3 mm, with blades that are slightly curved on the flat. Adventitia scissors are very sharp at the tip, with straight blades.

Dissection scissors are used for separating nerves and blood vessels from the surrounding tissue, without damaging them. When held closed, they are often also used as a dissecting probe. Dissection scissors are unsuitable for fine trimming work, since at high magnification the rounded tips prevent the surgeon from seeing exactly the tissue that is being cut, at the moment of cutting.

Adventitia scissors are used mostly at high magnification, for fine trimming and re-shaping of structures that have already been dissected free. They are also used for cutting fine micro suture material, 8/0 or less.

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SDA-type scissors

SDA-type scissors were specifically developed for fine dissection of the extremities and similar applications. Unlike other dissection scissors, the blades are **straight**. These scissors are held between thumb (top position) and index, stabilized with the middle and ring finger. You will immediately notice that handling the scissors in this position is easier than the usual "chuck grip". A further advantage is the better view of operating field.

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S&T®
Basic
informations
基本情報

Blade lengths

The blade length is measured from screw to the tip.

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S = short blades

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Westcott scissors

S&T® Westcott scissors feature finer blades than regular Westcotts, originally used in plastic surgery. The typical shape of the blades remains.

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Serrated scissors

The unique design of these handcrafted **S&T®** scissors ensures very smooth cutting action, even with the serrations.

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Vannas scissors

Originally designed for ophthalmic procedures, the vannas scissors are applied in vascular surgery, too. They feature finer, smaller and shorter blades.

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

Scissors by Blondeel

These scissors were designed by Professor Phillip Blondeel for use in breast reconstruction utilizing a perforator flap. While the ring handle configuration makes these suitable as preparation scissors for fine structures, the round handle, spring instrument configuration, with the curved blades, is extremely useful for fine trimming and dissection work. This instrument will perform tasks for which you might normally need different instruments.

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Special scissors

Flat handle

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Straight Dissection or curved Adventitia scissors? Our specially designed scissors meet all specific requirements!

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